
 **OMNIlife science, Inc.**  
480 Paramount Drive  
Raynham, MA 02767 USA  
www.omnils.com  
(508) 824-2444

#### LABEL SYMBOLS

	Use By
	Single Use Only
	See Instructions for Use
	Do Not Use if Package is Damaged
<b>Rx only</b>	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
<b>QTY</b>	Quantity
<b>STERILE R</b>	Sterilized Using Irradiation
<b>LOT</b>	Batch Code
<b>REF</b>	Catalogue Number
	Manufacturer

Medical Device Safety Service  
GmbH  
Schiffgraben 41  
30175 Hannover  
Germany

## OMNIBotics™ Tracker Kit

**Contents:**  
**Sterile Markers**  
**USB Drive**

#### DESCRIPTION

The accessories and instruments associated with the OMNIBotics applications make up a system designed for computer-assisted and robotic-assisted surgery. This packaging contains a set of sterile disposable markers and a non-sterile USB Drive. These markers are to be used with the instruments for the OMNIBotics computer-assisted and robotic-assisted surgery applications.

The appropriate number of markers (retro-reflective sensors) must be fixed onto the associated

OMNIBotics system instruments. The reflective surface of the markers allows a reference frame to be established, which can then be identified by an optical locator.

#### STERILITY AND PRODUCT HANDLING

Markers are provided sterile and intended for single use. They should always be stored unopened in their respective protective packaging. These sterile medical devices are delivered after sterilization by gamma ionizing radiation at a minimum dose of 25 kGy. These sterile devices must not be re-sterilized or reused.

Prior to use, inspect package for damage that may compromise sterility. If packaging has been opened or damaged, do not use and contact an OMNI representative. When removing the device from its packaging, relevant aseptic handling must be observed. Visually inspect each device prior to use for damage.

Sterile products, packaged in their original packaging, must be stored in a dry place away from light and at temperatures not exceeding 35°C (95°F). They must not be exposed to ionizing sources.

It is recommended that the surgeon count the markers after use.

#### MATERIAL

OMNIBotics System Markers are manufactured from HDPE polyethylene and 3M Scotchlite Reflective Tape.

#### INTENDED USE

The OMNIBotics System is intended for use during stereotaxic surgery to aid the surgeon in locating anatomical structures and aligning endoprostheses with the anatomical structures. A Surgical Technique can be obtained from an OMNI representative.

#### WARNINGS AND PRECAUTIONS

- OMNI devices must only be used by surgeons who have been trained in the surgical technique and are familiar with the instruments provided.

- The marker should be pushed firmly onto the clip until a clear “click” is both heard and felt, indicating that it has been correctly fixed to the rigid body. The marker must not be assembled or removed in the operating field, to avoid it falling into the incision.
- OMNI assumes no responsibility in the event of incorrect use of the instruments and accessories.
- The instruments must only be used with markers provided by OMNI. Otherwise, OMNI takes no responsibility for the surgical results.
- The surgeon must take care not to exert inappropriate stress on the device and must comply with the operating procedure described in the surgical technique.
- Prior to using the instrument system, the surgeon should give careful consideration to all aspects of the surgical intervention as well as the limitations of the implant and instruments.
- Sterile devices must not be re-sterilized. Resterilization of these devices will lead to deformation and imprecision in the localization system.
- The device is made of medical grade but not implant grade materials. It is recommended that the surgeon count the markers after use.
- Incorrect handling may render the instrument unsuitable for its intended use.
- The OMNIBotics System Markers that attach to the instruments are indicated for single-use only and must be detached and discarded according to the hospitals biohazard waste procedures.

All rights reserved. OMNIBotics™ is a trademark of OMNIlife science, Inc.  
Additional information may be obtained from OMNIlife science, Inc.

**SPANISH (ES)**

 **OMNI**life science, Inc.  
480 Paramount Drive  
Raynham, MA 02767 USA  
[www.omnils.com](http://www.omnils.com)  
(508) 824-2444

Medical Device Safety Service  
GmbH  
Schiffgraben 41  
30175 Hannover  
Germany

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