
OMNI life science, Inc.
 480 Paramount Dr.
 Raynham, MA 02767 USA
 www.omnils.com

(508)824-2444

European Authorized Representative:
 Medical Device Safety Service GmbH
 Schiffgraben 41
 30175 Hannover
 Germany

Cancellous Bone Screw

English- EN

SYMBOLS Glossary per ISO 15223-1

5.1.1		Medical Device Manufacturer
5.1.4		Use-By Date
5.4.2		Do not Re-use
5.4.4		See Instructions for Use
5.2.8		Do Not Use if Package is Damaged
	Rx only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
	QTY	Quantity
5.2.3		Sterilized Using Ethylene Oxide
5.1.5		Batch Code
5.1.6		Catalogue Number

PRODUCT HANDLING

Implants are provided sterile and should always be stored unopened in their respective protective containers. Prior to use, inspect package for damage which may compromise sterility. If packaging has been opened or damaged, contact manufacturer's representative. When unpacking the implant, verify the labeling for correct Cat. No. and size. When removing the implant from its packaging, the sterile technique must be observed. Protect the implant from contact with objects that may damage the surface finish. Inspect each implant prior to use for visual damage. Procedures for implanting and removal are available upon request.

DESCRIPTION

The Cancellous Bone Screws are intended to be used with the total joint systems provided by OMNIlife science. The Screws are provided in various lengths to accommodate the procedure. A surgical technique is available from OMNIlife science.

MATERIAL

- Titanium alloy (ASTM F 136)

INDICATIONS FOR USE

The Cancellous Bone Screw is intended to provide supplemental fixation when used with tibial and acetabular components in total knee and hip arthroplasty.

CONTRAINDICATIONS

Absolute contraindications include:

- Infection or sepsis or osteomyelitis;
- Insufficient bone structure or quality which may affect the stability of the implant;
- Rapid joint destruction or bone absorption;
- Skeletal immaturity;
- Muscular, ligamentous, neurological, vascular deficiencies or poor skin coverage, which may compromise the affected extremity;
- Alcoholism or other addictions;
- Sensitivity to the implant materials;
- High levels of physical activity (e.g. competitive sports, heavy physical labor);
- Obesity can produce loads on the prosthesis, which can lead to fixation failure or prosthesis breakage or fracture.

Relative contraindications include:

- Uncooperative patient or a patient with neurological disorders and incapable of following instruction;
- Metabolic disorders which may impair bone formation or bone quality;
- Distant foci of infections.

WARNINGS AND PRECAUTIONS

The Cancellous Bone Screws are provided sterile. Sterile technique should be used when opening the package.

Proper handling of the Screw before and during use is critical. Avoid bending, notching, or scratching the Screw as these may cause failure.

While these implants are not intended to withstand activity levels and loads of normal healthy bone, they are a means of restoring mobility and reducing pain for many patients.

In using total joint implants, the surgeon should be aware of the following:

- The correct selection of the modular implant components is extremely important. The potential for success in total joint replacement is increased by the selection of the proper size, shape and design of the implant. Total joint prostheses require careful seating and adequate bone support, and should be restricted to limited functional stress. The surgeon is to be thoroughly familiar with the implant, instruments, and surgical procedure prior to performing surgery.
- In selecting patients for total joint replacements, the following factors can be of extreme importance to the eventual success of the procedure:
 1. The patient's weight. An overweight or obese patient can produce loads on the prosthesis, which can lead to failure of the prosthesis. This becomes a major consideration when a small prosthesis must be used. Patients receiving knee joint replacements should be advised that the longevity of the implant may depend on their weight and level of activity.
 2. The patient's occupation or activity. If the patient is involved in an occupation or activity, which includes substantial walking, running, lifting, or muscle strain, the resultant forces can cause failure of the fixation, the device or both.

3. A condition of senility, mental illness or alcoholism. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions, leading to failure or other complications.
 4. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made prior to material selection or implantation.
 5. Certain degenerative diseases. In some cases, the progression of degenerative disease may be so advanced at the time of implantation that it may substantially decrease the expected useful life of the appliance. For such cases, total joint replacement can only be considered a delaying technique or temporary relief.
- The correct handling of the implant is extremely important. Care must be taken to protect surfaces from nicks and scratches which could become the focal point for failure. Do not tamper with the implant as contouring or bending of the implant may reduce its service life and may cause immediate or eventual failure under load.
 - A surgical implant should not be reused. Even though a used implant may appear undamaged, it may have small defects and internal stress patterns, which may lead to failure. Use only new prosthesis of the current design.
 - Resterilization of the device is not recommended.
 - Bone excision should be limited to the amount necessary to accommodate the implants. Prior to closure, the surgical site should be thoroughly cleaned of bone chips, bone cement or other detritus that may cause a third body wear problem. Range of motion should be checked for impingement or instability.
 - Postoperative care is important. The patient should be instructed on the limitations of these devices and should be cautioned regarding load-bearing, ranges of motion, and activity levels permissible. Excessive physical activity and trauma affecting the replaced joint have been implicated in premature failure by loosening, fracture and/or wear of the prosthesis implant. Early load-bearing should be carefully controlled. The patient should be advised to report any related pain, decrease in range of motion, swelling, fever, and unusual incidences.
 - This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

POSSIBLE ADVERSE EFFECTS

The possible adverse effects are similar to those occurring with any hip or knee replacement and include the following:

- Loosening or migration of components due to trauma and/or loss of fixation.
- Histiocytic granuloma formation and osteolysis around the implant due to wear debris.
- Fracture of the implant as the result of strenuous activity, improper alignment, inadequate fixation or extreme duration of service.
- Urological complications, especially urinary retention and infection.
- Other complications associated with general surgery, drugs or ancillary devices used, blood, etc.

Intraoperative and early postoperative complications can include:

- Damage to blood vessels;
- Temporary or permanent neuropathies;
- Traumatic arthrosis of the knee from Intraoperative positioning of the extremity;
- Cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction;
- Hematoma;
- Delayed wound healing;
- Infection
- Femoral Perforation

Late postoperative complications can include:

- Aggravated problems of the hip or ankle of the affected limb or contralateral extremity by leg length discrepancy or muscle deficiency;
- Periarticular calcification or ossification, with or without impediment to joint mobility;
- Infection;
- Trochanteric non-union due to inadequate reattachment and/or early weight bearing.

CAUTION

Disposal of implants should be carried out using the hospital's standard method for non-biodegradable non-combustible medical waste.

MRI SAFETY INFORMATION

The implants have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the implants in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

All rights reserved. OMNI Knee™ and Apex Knee™ System are trademarks of OMNIlife science, Inc.

Additional information about the Apex Knee™ System may be obtained from OMNI life science, Inc.

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