

Manufactured by:


OMNIlife Science, Inc. (OMNI)

480 Paramount Drive
Raynham, MA 02767 USA
www.omnils.com

For technical support call:

☎ : 508-824-2444

📞 : 508-822-6030



Medical Device Safety Service GmbH
Schiffgraben 41
30175 Hannover
Germany

INTENDED USE

The iBlock Sawing Guide Very Short is intended to function exactly as the existing iBlock Sawing Guides (#4146-5200 and #4146-5300) while offering a smaller device footprint. This smaller sawing guide can be used in the case of a smaller incision, where the larger sawing guides might contact soft tissues that will impede on the devices mobility.

REFERENCES

- IFU-036: OMNI ART KNEE APPLICATION AND IBLOCK (OMNIBotics Total Knee) IFU
- IFU-039: OMNI ART AND IBLOCK SYSTEM INSTRUMENTATION IFU

INSTRUCTIONS FOR USE

The Very Short iBlock Sawing Guide's instructions for use can directly mirror section 4.16 of IFU-036 with some minor changes:

- In all places where the IFU reads "iBlock Cutting Guide (#4146-5200 or #4146-5300)", it should instead read "iBlock Cutting Guide (#4146-5200 or #4146-5300 or #4146-5400)".

